

## Cell Culture Flasks

### I. APPLICATIONS

1. Suitable for cell and tissue culture

### II. FEATURES

1. Constructed of high quality medical grade polystyrene
2. Sterile
3. Non-pyrogenic
4. Non-cytotoxic
5. DNase, RNase, and human DNA free
6. Stackable
7. Surface Treatment

### III. DESIGN

1. Cell culture surface area of 25 cm<sup>2</sup>, 75 cm<sup>2</sup>, or 175 cm<sup>2</sup>
2. Plug cap and filter cap options available for all flask sizes
3. Short, wide, angled neck design for easy access

### IV. PACKAGING

1. Bulk, 5 Pieces/bag, 200 Pieces/case. (07-8025, 07-9025)
2. Bulk, 5 Pieces/bag, 100 Pieces/case. (07-8075, 07-9075)
3. Bulk, 5 Pieces/bag, 40 Pieces/case. (07-8175, 07-9175)

### V. DESCRIPTION

Item	Cap type	Culture Area (cm <sup>2</sup> )	Working Volume (ml)	Total Volume (ml)	Packaging
07-8025	Filter	25	7	60	5 Pieces/bag, 200 Pieces/case
07-8075	Filter	75	25	250	5 Pieces/bag, 100 Pieces/case
07-8175	Filter	175	50	650	5 Pieces/bag, 40 Pieces/case
07-9025	plugged	25	7	60	5 Pieces/bag, 200 Pieces/case
07-9075	plugged	75	25	250	5 Pieces/bag, 100 Pieces/case
07-9175	plugged	175	50	650	5 Pieces/bag, 40 Pieces/case

### VI. STERILITY

1. Sterile

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard

**VIII. RESTRICTIONS**

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

**IX. PRODUCT PARAMETER**

Dimensions=mm

Cat #	Description	A	B	C	D	E
<b>07-8025</b>	Cell Culture Flask,Filter Cap	50.10	72.00	92.38	14.50	24.60
<b>07-8075</b>	Cell Culture Flask,Filter Cap	84.20	124.10	151.04	19.00	32.00
<b>07-8175</b>	Cell Culture Flask,Filter Cap	116.60	199.08	226.69	20.50	34.72
<b>07-9025</b>	Cell Culture Flask,Plug Cap	50.10	72.00	92.38	14.50	24.60
<b>07-9075</b>	Cell Culture Flask,Plug Cap	84.20	124.10	151.04	19.00	32.00
<b>07-9175</b>	Cell Culture Flask,Plug Cap	116.60	199.08	226.69	20.50	34.72

