

Gloves

I. APPLICATIONS

1. Suitable for all types of medical examination, home care and diagnostic and laboratory use

II. GLOVES FEATURES

Powder Free Latex Disposable Glove Features

1. Powder-free, single chlorinated Latex Gloves to provide good dry hand donning capability
2. Full textured surface enhances grip for efficient instrument handling
3. Good grip to handle dry instruments when moderate amount of handling is required
4. Good Thickness Distribution for tactile sensitivity
5. Suitable for all types of medical examination, home care and diagnostic and laboratory use
6. Natural to Light Yellow Color

Powder Free Nitrile Disposable Glove Features

1. Non-added components (sulphur and accelerator chemicals) to minimize exposure to Type IV hypersensitivity
2. Non-latex formulation eliminates exposure to Type 1 latex protein sensitivity
3. Nitrile material promotes better puncture and chemical/ chemo resistance properties
4. Low Dermatitis Potential Claim
5. Finger textured with improved working grip to enable efficient instrument handling in wet and dry condition
6. Beaded cuff design to don gloves efficiently and securely
7. Extra length
8. Dark Blue Color

III. GLOVES DESIGN

1. Two kinds of gloves: Latex and Nitrile
2. Non-sterile
3. Ambidextrous
4. Biocompatibility Tested (ISO 10993)
5. Patent Pending Formulation

IV. PACKAGING

100 Pieces/pack, 10 Packs/case

V. DESCRIPTION

Catalog #	Description	Size	Color	Sterile	Unit
97-1111	Latex	XS	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1112	Latex	S	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1113	Latex	M	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1114	Latex	L	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1115	Latex	XL	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-6111	Nitrile	XS	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6112	Nitrile	S	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6113	Nitrile	M	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6114	Nitrile	L	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6115	Nitrile	XL	Dark Blue	No	100 Pieces/pack, 10 Packs/case

VI. STERILITY

- No

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Catalog #	Size	A
97-1111	XS	70
97-1112	S	80
97-1113	M	90
97-1114	L	100
97-1115	XL	110
97-6111	XS	70
97-6112	S	80
97-6113	M	90
97-6114	L	100
97-6115	XL	110

